

Anti-Lig1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P18858

Item Number: CM0022692



Introduction

Rabbit polyclonal antibody against human DNA ligase 1 (LIG1). Validated for WB, IF/ICC, ELISA, with cross-reactivity to human and mouse. LIG1 plays a key role in DNA repair, replication, and recombination.

[Learn More](#)

Parameter	Value
Product Name	LIG1 Rabbit pAb
Remarks/Alias	LIG1; IMD96; hLig1; LIG1
Species	Human
GeneID (Human)	3978
GeneID	3978
Immunogen	Recombinant fusion protein containing amino acids 565-790 of human LIG1 (NP_000225.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P18858
Protein Weight	102kDa
Shipping	Ice bag