

## Anti-Lilra2/Cd85H Rabbit Monoclonal Antibody For If/Icc, Fc, Elisa - Q8N149

Item Number: CM0009268



### Introduction

LILRA2/CD85h rabbit monoclonal antibody validated for IF/ICC, FC, ELISA. Reacts specifically with human LILRA2 (UniProt Q8N149). Suitable for innate immunity and microbial infection research applications.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | LILRA2/CD85h Rabbit mAb                                                                                                     |
| Remarks/Alias    | ILT1; LIR7; CD85H; LIR-7                                                                                                    |
| Species          | Human                                                                                                                       |
| GeneID (human)   | 11027                                                                                                                       |
| GeneID           | 11027                                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 24-449 of human LILRA2/CD85h (NP_001124389.2) |
| Source           | Rabbit                                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                                       |
| Application      | IF/ICC, FC, ELISA                                                                                                           |
| Cross Reactivity | Human                                                                                                                       |
| SWISS            | Q8N149                                                                                                                      |
| Protein Weight   | 53kDa/51kDa/48kDa/47kDa                                                                                                     |
| Shipping         | Ice bag                                                                                                                     |