


Anti-Lilrb1/Cd85j/Ilt2 Monoclonal Antibody For Wb, Fc, Elisa - Q8Nhl6

Item Number: CM0006854



Introduction

Rabbit monoclonal antibody against human LILRB1 (CD85j/ILT2) validated for Western blot, flow cytometry, and ELISA. Inhibitory receptor for MHC class I; suitable for immune checkpoint and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	LILRB1/CD85j/ILT2 Rabbit mAb
Remarks/Alias	ILT2; LIR1; MIR7; PIRB; CD85j; ILT-2; LIR-1; MIR-7; PIR-B; LILRB1/CD85j/ILT2
Species	Human
GeneID	10859
Immunogen	Recombinant protein of human LILRB1/CD85j/ILT2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	Q8NHL6
Protein Weight	49KDa/70KDa/71KDa
Shipping	Ice bag