

Anti-Lilrb2/Cd85D/Ilt4 Rabbit Polyclonal Antibody For Wb/Elisa - Q8N423

Item Number: CM0013087



Introduction

Rabbit polyclonal antibody against human LILRB2 (CD85d/ILT-4), validated for WB and ELISA. Ideal for studying immune tolerance, HLA-G interaction, and receptor signaling in monocytes and dendritic cells.

[Learn More](#)

Parameter	Value
Product Name	LILRB2/CD85d/ILT4 Rabbit pAb
Remarks/Alias	ILT4; LIR2; CD85D; ILT-4; LIR-2; MIR10; MIR-10; LILRB2/CD85d/ILT4
Species	Human
GeneID (Human)	10288
GeneID	10288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 479-598 of human LILRB2/CD85d/ILT4 (NP_005865.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N423
Protein Weight	65kDa
Shipping	Ice bag