


Anti-Lilrb3 Polyclonal Antibody For Wb, Elisa - 075022

Item Number: CM0020716



Introduction

Rabbit polyclonal antibody against human LILRB3 (ILT5/CD85a), validated in WB and ELISA. Cross-reacts with mouse. LILRB3 is an inhibitory receptor on B cells and monocytes, down-regulating activation. Suitable for immune research.

[Learn More](#)

Parameter	Value
Product Name	LILRB3 Rabbit pAb
Remarks/Alias	HL9; ILT5; LIR3; PIRB; CD85A; ILT-5; LIR-3; PIR-B; LILRA6; LILRB3
Species	Human
Gene ID (Human)	11025
Gene ID	11025
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075022
Protein Weight	69kDa
Shipping	Ice bag