

Anti-Lilrb4 Polyclonal Antibody For Wb, Ihc, If, Elisa - Q8Nhj6

Item Number: CM0025419



Introduction

Rabbit polyclonal antibody against human LILRB4 (ILT3/CD85K), suitable for WB, IHC-P, IF-P, ELISA; detects human, mouse, rat. Targets an inhibitory immune receptor for tolerance and immuno-oncology studies.

[Learn More](#)

Parameter	Value
Product Name	LILRB4 Rabbit pAb
Remarks/Alias	B4; ILT3; LIR5; CD85K; ILT-3; LIR-5; LILRB4
Species	Human
Gene ID	11006
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NHJ6
Protein Weight	49kDa
Shipping	Ice bag