

Anti-Lim Kinase 1 (Limk1) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P53667

Item Number: CM0006715



Introduction

IVD-grade rabbit monoclonal antibody against LIM kinase 1 (LIMK1) for Western Blot, IHC-P, and ELISA. Detects human, mouse, and rat LIMK1. Ideal for actin cytoskeleton and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	LIM kinase 1(LIMK1) Rabbit mAb
Remarks/Alias	LIMK1; LIMK; LIMK-1; LIM domain kinase 1; LIM kinase 1(LIMK1)
Species	Human
Gene ID (Human)	3984
Gene ID	3984
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 630-647 of human LIM kinase 1(LIMK1) (NP_002305.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53667
Protein Weight	33kDa/68kDa/70kDa/72kDa
Shipping	Ice bag