

Anti-Lim2 Rabbit Polyclonal Antibody For Wb And Elisa - P55344

Item Number: CM0024978



Introduction

High-quality rabbit polyclonal antibody targeting human LIM2 (MP19), validated for WB and ELISA. Suitable for lens development and cataractogenesis research in human, mouse, and rat samples.

[Learn More](#)

Parameter	Specification
Product Name	LIM2 Rabbit pAb
Remarks/Alias	MP17, MP19, CTRCT19, LIM2
Species	Human
GeneID (Human)	3982
GeneID	3982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-110 of human LIM2 (NP_085915.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55344
Protein Weight	20kDa
Shipping	Ice bag