

## Anti-Lima1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uhb6

Item Number: CM0017022



### Introduction

High-quality rabbit polyclonal antibody specific for human LIMA1 (EPLIN), validated in WB, IHC-P, and ELISA. Detects mouse LIMA1. Supports research on actin dynamics, primary cilium, and cholesterol absorption.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | LIMA1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | EPLIN; LDLCQ8; SREBP3; LIMA1                                                                                       |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 51474                                                                                                              |
| GeneID           | 51474                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 560-759 of human LIMA1 (NP_057441.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse                                                                                                       |
| SWISS            | Q9UHB6                                                                                                             |
| Protein Weight   | 85kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |