

Anti-Limch1 Polyclonal Antibody For Wb, Elisa - Q9Upq0

Item Number: CM0022121

Introduction

High-quality anti-LIMCH1 rabbit polyclonal antibody for Western blot and ELISA applications. Validated in human samples. Targets actin stress fiber associated protein regulating cell migration and adhesion.

[Learn More](#)

Parameter	Value
Product Name	LIMCH1 Rabbit pAb
Remarks/Alias	LMO7B; LIMCH1A; LIMCH1
Species	Human
GeneID (Human)	22998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-450 of human LIMCH1 (NP_001106189.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UPQ0
Protein Weight	122kDa
Shipping	Ice bag