


Anti-Lims1 Polyclonal Antibody For Wb, If/Icc, Elisa - P48059

Item Number: CM0017114



Introduction

Rabbit polyclonal antibody against human LIMS1 (PINCH1), validated for WB, IF/ICC, ELISA. Suitable for human and mouse samples. Ideal for cell adhesion and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	LIMS1 Rabbit pAb
Remarks/Alias	PINCH; PINCH1; PINCH-1; LIMS1
Species	Human
GeneID (Human)	3987
GeneID	3987
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human LIMS1 (NP_001180417.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48059
Protein Weight	37kDa
Shipping	Ice bag