

Anti-Lims1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P48059

Item Number: CM0030815



Introduction

High-quality rabbit polyclonal antibody targeting human LIMS1 (PINCH-1). Validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for focal adhesion and cell migration studies. 37kDa.

[Learn More](#)

Parameter	Value
Product Name	LIMS1 Rabbit pAb
Remarks/Alias	PINCH; PINCH1; PINCH-1; LIMS1
Species	Human
Gene ID (Human)	3987
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human LIMS1 (NP_001180417.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48059
Protein Weight	37kDa
Shipping	Ice bag