

Anti-Lin28A Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H9Z2

Item Number: CM0007221



Introduction

Recombinant rabbit monoclonal antibody targeting human LIN28A, validated for Western blot, IHC-P, and ELISA. Detects a 23 kDa protein and cross-reacts with mouse and rat. Ideal for pluripotency and cancer research.

[Learn More](#)

Parameter	Value
Product Name	LIN28A Rabbit PolymAb®
Remarks/Alias	CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A
Species	Human
GeneID (Human)	79727
GeneID	79727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human LIN28A. (NP_078950.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9Z2
Protein Weight	23kDa
Shipping	Ice bag