

Anti-Lin28A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H9Z2

Item Number: CM0027786



Introduction

High-quality anti-LIN28A rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Targets human and mouse LIN28A, a key regulator of pluripotency and miRNA biogenesis. Ideal for stem cell and oncology research.

[Learn More](#)

Parameter	Value
Product Name	LIN28A Rabbit pAb
Remarks/Alias	CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A; LIN28A
Species	Human
GeneID (Human)	79727
GeneID	79727
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-209 of human LIN28A (NP_078950.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H9Z2
Protein Weight	23kDa
Shipping	Ice bag