


Anti-Lin28B Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q6Zn17

Item Number: CM0012596



Introduction

Rabbit monoclonal antibody targeting human LIN28B, validated for WB, IP, IHC-P, and ELISA. LIN28B is a key regulator of microRNA biogenesis, crucial for stem cell pluripotency and cancer progression. Ideal for diagnostic and research use.

[Learn More](#)

Parameter	Value
Product Name	LIN28B Rabbit mAb
Remarks/Alias	CSDD2; LIN28B
Species	Human
Gene ID (Human)	389421
Gene ID	389421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human LIN28B (Q6ZN17)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q6ZN17
Protein Weight	27kDa
Shipping	Ice bag