

Anti-Lin7C Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nup9

Item Number: CM0023677



Introduction

Anti-LIN7C rabbit polyclonal antibody validated for Western blotting and ELISA. Reacts with human, mouse, and rat LIN7C (Q9NUP9). Suitable for studying cell polarity, membrane trafficking, and synaptic organization.

[Learn More](#)

Parameter	Value
Product Name	LIN7C Rabbit pAb
Remarks/Alias	MALS3; VELI3; LIN-7C; MALS-3; LIN-7-C; LIN7C
Species	Human
GeneID (Human)	55327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human LIN7C (NP_060832.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NUP9
Protein Weight	22kDa
Shipping	Ice bag