

Anti-Lin9 Polyclonal Antibody For Wb, Elisa - Q5Tka1

Item Number: CM0024434



Introduction

Rabbit polyclonal antibody targeting human LIN9 (TGS, BARA), validated for WB and ELISA. Suitable for human and mouse samples. Ideal for cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	LIN9 Rabbit pAb
Remarks/Alias	TGS; BARA; TGS1; TGS2; Lin-9; BARPsv; LIN9
Species	Human
GeneID (Human)	286826
GeneID	286826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-470 of human LIN9 (NP_775106.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5TKA1
Protein Weight	62kDa
Shipping	Ice bag