

Anti-Lingo1 Polyclonal Antibody For Wb, Elisa - Q96Fe5

Item Number: CM0024874



Introduction

LINGO1 Rabbit pAb is a rabbit polyclonal antibody targeting human LINGO1 (Q96FE5), validated for Western blot and ELISA, with cross-reactivity to mouse, suitable for studies of axonal regeneration and myelination.

[Learn More](#)

Parameter	Value
Product Name	LINGO1 Rabbit pAb
Remarks/Alias	LERN1, MRT64, LRRN6A, UNQ201, LINGO1
Species	Human
GeneID (Human)	84894
GeneID	84894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-240 of human LINGO1 (NP_116197.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96FE5
Protein Weight	70kDa
Shipping	Ice bag