

Anti-Lins1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Ng48

Item Number: CM0014439

Introduction

Anti-LINS1 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, and rat LINS1 protein (86 kDa). Ideal for Wnt signaling and intellectual disability research.

[Learn More](#)

Parameter	Value
Product Name	LINS1 Rabbit pAb
Remarks/Alias	LINS; MRT27; WINS1; LINS1
Species	Human
GeneID	55180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 541-700 of human LINS1 (NP_001035706.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NG48
Protein Weight	86kDa
Shipping	Ice bag