

Anti-Lipc (Hepatic Triacylglycerol Lipase) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P11150

Item Number: CM0017075

Introduction

High-quality rabbit polyclonal antibody targeting human LIPC (hepatic triacylglycerol lipase). Validated for WB, IF/ICC, and ELISA. Ideal for lipid metabolism, HDL, and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	LIPC Rabbit pAb
Remarks/Alias	HL; HTGL; LIPH; HDLCQ12; LIPC
Species	Human
Gene ID (Human)	3990
Gene ID	3990
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-499 of human LIPC (NP_000227.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P11150
Protein Weight	56kDa
Shipping	Ice bag