


Anti-Lipf Rabbit Polyclonal Antibody For Wb And Elisa - P07098

Item Number: CM0030045



Introduction

LIPF rabbit polyclonal antibody for WB, ELISA applications; recognizes human, mouse, rat. Targets gastric triacylglycerol lipase (P07098); secreted enzyme for lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	LIPF Rabbit pAb
Remarks/Alias	GL; HGL; HLAL; LIPF
Species	Human
Gene ID (Human)	8513
Gene ID	8513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 219-398 of human LIPF (NP_004181.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07098
Protein Weight	45kDa
Shipping	Ice bag