

Anti-Lipg Polyclonal Antibody For Wb, Elisa - Q9Y5X9

Item Number: CM0023184

Introduction

Rabbit polyclonal antibody against human Endothelial Lipase (LIPG). Validated for WB and ELISA; cross-reacts with mouse. Ideal for lipid metabolism and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	LIPG Rabbit pAb
Remarks/Alias	EL; EDL; PRO719; LIPG
Species	Human
Gene ID	9388
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 251-500 of human LIPG (NP_006024.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5X9
Protein Weight	57kDa
Shipping	Ice bag