

Anti-Lipin 1 Rabbit Polyclonal Antibody For Wb And Elisa - Q14693

Item Number: CM0016719

Introduction

Rabbit polyclonal antibody targeting human Lipin-1 (LPIN1), validated for Western blot and ELISA. Cross-reacts with mouse and rat, ideal for lipid metabolism and adipogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Lipin 1 Rabbit pAb
Remarks/Alias	PAP1; Lipin 1
Species	Human
GeneID	23175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 261-414 of human Lipin 1 (NP_663731.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14693
Protein Weight	99kDa
Shipping	Ice bag