

Anti-Lipin 1 Rabbit Polyclonal Antibody For Wb, Elisa - Q14693

Item Number: CM0026227

Introduction

Rabbit polyclonal antibody targeting human Lipin 1 (LPIN1), validated for Western blot and ELISA. Detects human, mouse and rat species. Useful for studying lipid metabolism, adipocyte differentiation, and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	Lipin 1 Rabbit pAb
Remarks/Alias	PAP1; Lipin 1
Species	Human
GeneID (human)	23175
GeneID	23175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 261-414 of human Lipin 1 (NP_663731.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14693
Protein Weight	99kDa
Shipping	Ice bag