


Anti-Lipocalin-2/NGal Rabbit Polyclonal Antibody For Wb, Elisa - P80188

Item Number: CM0027369



Introduction

Rabbit polyclonal antibody against human Lipocalin-2/NGAL (SWISS P80188). Validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for iron trafficking, immunity, and biomarker research.

[Learn More](#)

Parameter	Specification
Product Name	Lipocalin-2/NGAL Rabbit pAb
Remarks/Alias	p25; 24p3; MSFI; NGAL; Lipocalin-2/NGAL
Species	Human
GeneID (human)	3934
GeneID	3934
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P80188
Protein Weight	23kDa
Shipping	Ice bag