

Anti-Lipoprotein Lipase (Lpl) Rabbit Monoclonal Antibody For Wb, Elisa - P06858

Item Number: CM0008438



Introduction

Recombinant rabbit monoclonal antibody against human Lipoprotein lipase (LPL, P06858). Validated for WB and ELISA, cross-reacts with mouse. Suitable for lipid metabolism and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	Lipoprotein lipase (LPL) Rabbit mAb
Remarks/Alias	LIPD; HDLCQ11; Lipoprotein lipase (LPL)
Species	Human
GeneID (Human)	4023
GeneID	4023
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 376-475 of human Lipoprotein lipase (LPL) (P06858).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06858
Protein Weight	53kDa
Shipping	Ice bag