

Anti-Lipoprotein Lipase (Lpl) Rabbit Polyclonal Antibody For Wb, Elisa - P06858

Item Number: CM0022265



Introduction

Rabbit polyclonal antibody against human Lipoprotein lipase (LPL), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for lipid metabolism and cardiovascular research. SWISS: P06858.

[Learn More](#)

Parameter	Value
Product Name	Lipoprotein lipase (LPL) Rabbit pAb
Remarks/Alias	LIPD; HDLCQ11; Lipoprotein lipase (LPL)
Species	Human
Gene ID (Human)	4023
Gene ID	4023
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 124-475 of human Lipoprotein lipase (Lipoprotein lipase (LPL)) (NP_000228.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06858
Protein Weight	53kDa
Shipping	Ice bag