


Anti-Lis1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P43034

Item Number: CM0011212



Introduction

Anti-LIS1 (PAFAH1B1) rabbit monoclonal antibody for WB, IF/ICC, IP, ELISA. Validated in human, mouse, rat; detects 47kDa protein involved in dynein regulation and neuronal migration.

[Learn More](#)

Parameter	Value
Product Name	LIS1 Rabbit mAb
Remarks/Alias	MDS; LIS1; LIS2; MDCR; NudF; PAFAH
Species	Human
Gene ID (Human)	5048
Gene ID	5048
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LIS1 (P43034).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43034
Protein Weight	47kDa
Shipping	Ice bag