


Anti-Lkb1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q15831

Item Number: CM0021725



Introduction

LKB1 rabbit polyclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Targets tumor suppressor STK11 (Q15831), ideal for metabolism, polarity, and DNA damage studies.

[Learn More](#)

Parameter	Value
Product Name	LKB1 Rabbit pAb
Remarks/Alias	PJS; LKB1; hLKB1
Species	Human
GeneID (Human)	6794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 134-433 of human LKB1 (NP_000446.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15831
Protein Weight	49 kDa
Shipping	Ice bag