


Anti-Lkb1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q15831

Item Number: CM0004747



Introduction

Human LKB1 (STK11) rabbit monoclonal antibody validated for Western blot, IHC-P, IP, and ELISA. Detects endogenous LKB1 across human, mouse, and rat samples. Ideal for metabolism, cell polarity, and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] LKB1 Rabbit mAb
Remarks/Alias	PJS; LKB1; hLKB1; B1
Species	Human
Gene ID (Human)	6794
Gene ID	6794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 134-433 of human LKB1. (NP_000446.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15831
Protein Weight	49kDa
Shipping	Ice bag