

Anti-LLgl1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15334

Item Number: CM0018751



Introduction

Anti-LLGL1 rabbit polyclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat LLGL1, a 115 kDa cell polarity protein involved in tumor suppression and neurodevelopment.

[Learn More](#)

Parameter	Value
Product Name	LLGL1 Rabbit pAb
Remarks/Alias	DLG4; HUGL; LLGL; Lgl1; Mgl1; HUGL1; HUGL-1; LLGL1
Species	Human
GeneID (Human)	3996
GeneID	3996
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human LLGL1 (NP_004131.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15334
Protein Weight	115kDa
Shipping	Ice bag