


Anti-LLph Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Brt6

Item Number: CM0023814



Introduction

LLPH Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human and rat. Suitable for studies of dendritic and spine growth and synaptic transmission.

[Learn More](#)

Parameter	Value
Product Name	LLPH Rabbit pAb
Remarks/Alias	hLLP; C12orf31; LLPH
Species	Human
Gene ID (Human)	84298
Gene ID	84298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-129 of human LLPH (NP_115714.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BRT6
Protein Weight	15kDa
Shipping	Ice bag