


Anti-Lman1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P49257

Item Number: CM0011186



Introduction

LMAN1 Rabbit Monoclonal Antibody validated for WB, IHC-P, and ELISA; detects human, mouse, and rat ERGIC-53, a cargo receptor for ER-to-Golgi transport. Research-use-only.

[Learn More](#)

Parameter	Specification
Product Name	LMAN1 Rabbit mAb
Remarks/Alias	MR60, gp58, F5F8D, FMFD1, MCFD1, ERGIC53, ERGIC-53, LMAN1
Species	Human
Gene ID (Human)	3998
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LMAN1 (P49257).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49257
Protein Weight	54kDa
Shipping	Ice bag