


Anti-Lman1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49257

Item Number: CM0016531



Introduction

Rabbit polyclonal antibody against human LMAN1 (ERGIC-53, MR60) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Supports research on intracellular lectins and ER-Golgi cargo transport.

[Learn More](#)

Parameter	Value
Product Name	LMAN1 Rabbit pAb
Remarks/Alias	MR60; gp58; F5F8D; FMFD1; MCFD1; ERGIC53; ERGIC-53; LMAN1
Species	Human
GeneID (Human)	3998
GeneID	3998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-480 of human LMAN1 (NP_005561.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49257
Protein Weight	54kDa
Shipping	Ice bag