

## Anti-Lman2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q12907

Item Number: CM0016507

### Introduction

Rabbit polyclonal antibody against human LMAN2 (VIP36), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. LMAN2 is an intracellular lectin involved in glycoprotein trafficking in the early secretory pathway. Ideal for cell biology research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | LMAN2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | GP36B; VIP36; C5orf8; LMAN2                                                                                        |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 10960                                                                                                              |
| GeneID           | 10960                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 190-320 of human LMAN2 (NP_006807.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, ELISA                                                                                                  |
| Cross Reactivity | Mouse                                                                                                              |
| SWISS            | Q12907                                                                                                             |
| Protein Weight   | 40kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |