

Anti-Lman2L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H0V9

Item Number: CM0019716

Introduction

Rabbit polyclonal antibody targeting human LMAN2L (VIPL), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for investigating ER-to-Golgi glycoprotein trafficking and ERGIC-53 regulation.

[Learn More](#)

Parameter	Value
Product Name	LMAN2L Rabbit pAb
Remarks/Alias	VIPL; MRD69; MRT52; LMAN2L
Species	Human
Gene ID (Human)	81562
Gene ID	81562
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 189-348 of human LMAN2L (NP_110432.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0V9
Protein Weight	40kDa
Shipping	Ice bag