

Anti-Lmbrd1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nun5

Item Number: CM0013409

Introduction

Rabbit polyclonal antibody targeting human LMBRD1, validated for Western blot, IF/ICC, and ELISA. Detects human and mouse LMBRD1, useful for studying cobalamin transport and insulin receptor trafficking.

[Learn More](#)

Parameter	Value
Product Name	LMBRD1 Rabbit pAb
Remarks/Alias	NESI; LMBD1; MAHCF; C6orf209; LMBRD1
Species	Human
GeneID (Human)	55788
GeneID	55788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human LMBRD1 (NP_060838.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NUN5
Protein Weight	61kDa
Shipping	Ice bag