

Anti-Lmo1 (Rhombotin-1) Rabbit Polyclonal Antibody For Wb, Elisa - P25800

Item Number: CM0025953

Introduction

Rabbit polyclonal antibody targeting human LMO1 (Rhombotin-1), validated for WB and ELISA. Suitable for studies on neural gene regulation, LIM domain proteins, and T-cell leukemia. Specific to human, with an ~18 kDa predicted band.

[Learn More](#)

Parameter	Value
Product Name	LMO1 Rabbit pAb
Remarks/Alias	TTG1; RBTN1; RHOM1; LMO1
Species	Human
GeneID (Human)	4004
GeneID	4004
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human LMO1 (NP_002306.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P25800
Protein Weight	18kDa
Shipping	Ice bag