

Anti-Lmo2 Monoclonal Antibody For Wb, Elisa - P25791

Item Number: CM0004141



Introduction

High-quality Anti-LMO2 Rabbit Monoclonal Antibody for WB and ELISA. Targets human LMO2 (18 kDa), a nuclear transcription regulator essential for erythroid development. Ideal for hematopoiesis and leukemia research from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	LMO2 Rabbit mAb
Remarks/Alias	TTG2; LMO-2; RBTN2; RHOM2; RBTN1; LMO2
Species	Human
GeneID (Human)	4005
GeneID	4005
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P25791
Protein Weight	18kDa
Shipping	Ice bag