


Anti-Lmo4 Rabbit Monoclonal Antibody For Wb And Elisa - P61968

Item Number: CM0007859



Introduction

Rabbit monoclonal antibody targeting human LMO4 (LIM domain transcription factor LMO4), suitable for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying transcription regulation and motor neuron development.

[Learn More](#)

Parameter	Value
Product Name	LMO4 Rabbit mAb
Remarks/Alias	LMO4
Species	Human
Gene ID	8543
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 66-165 of human LMO4 (P61968)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61968
Protein Weight	18kDa
Shipping	Ice bag