


Anti-Lmo4 Rabbit Polyclonal Antibody For If/Icc And Elisa - P61968

Item Number: CM0027672



Introduction

Rabbit anti-LMO4 polyclonal antibody validated for IF/ICC and ELISA. Reacts with human, mouse, rat LMO4. Ideal for transcription regulation, motor neuron identity, and oncoprotein studies. SWISS: P61968.

[Learn More](#)

Parameter	Value
Product Name	LMO4 Rabbit pAb
Remarks/Alias	LMO4
Species	Human
Gene ID (Human)	8543
Gene ID	8543
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human LMO4 (NP_006760.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61968
Protein Weight	18kDa
Shipping	Ice bag