


Anti-Lmod1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P29536

Item Number: CM0030751



Introduction

Rabbit polyclonal antibody against human LMOD1, validated in WB, IHC-P, ELISA. Targets leiomodins-1, an actin nucleation factor critical for smooth muscle contractility. Suitable for research in muscle biology and autoimmunity.

[Learn More](#)

Parameter	Specification
Product Name	LMOD1 Rabbit pAb
Remarks/Alias	1D; D1; 64kD; MMIHS3; SMLMOD; SM-LMOD; LMOD1
Species	Human
GeneID (Human)	25802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-187 of human LMOD1 (NP_036266.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P29536
Protein Weight	67kDa
Shipping	Ice bag