

Anti-Lmx1B Polyclonal Antibody For Wb, Elisa - O60663

Item Number: CM0016937

Introduction

Rabbit polyclonal antibody against human LMX1B, validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for studying podocyte transcription regulation and limb development.

[Learn More](#)

Parameter	Value
Product Name	LMX1B Rabbit pAb
Remarks/Alias	NPS1; FSGS10; LMX1.2; LMX1B
Species	Human
GeneID (Human)	4010
GeneID	4010
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human LMX1B (NP_001167618.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60663
Protein Weight	45kDa
Shipping	Ice bag