

## Anti-Lmx1B Rabbit Polyclonal Antibody For Wb, Elisa - O60663

Item Number: CM0027315

### Introduction

Recombinant-based rabbit polyclonal antibody targeting human LMX1B transcription factor. Validated for WB and ELISA, cross-reacts with rat. Suitable for podocyte biology and limb development research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | LMX1B Rabbit pAb                                                                                                   |
| Remarks/Alias    | NPS1; FSGS10; LMX1.2; LMX1B                                                                                        |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 4010                                                                                                               |
| GeneID           | 4010                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 136-395 of human LMX1B (NP_002307.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Rat                                                                                                                |
| SWISS            | O60663                                                                                                             |
| Protein Weight   | 45kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |