

Anti-Lnx1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tbb1

Item Number: CM0027058



Introduction

LN1X1 rabbit polyclonal antibody suitable for Western blot (WB) and ELISA; recognizes human, mouse, and rat LN1X1, an E3 ubiquitin-protein ligase that mediates NUMB degradation and endocytic scaffolding.

[Learn More](#)

Parameter	Value
Product Name	LN1X1 Rabbit pAb
Remarks/Alias	LN1X; MPDZ; PDZRN2; LN1X1
Species	Human
GeneID (human)	84708
GeneID	84708
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 501-594 of human LN1X1 (NP_001119800.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBB1
Protein Weight	81kDa
Shipping	Ice bag