

Anti-Lonp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Wa8

Item Number: CM0022885

Introduction

Rabbit polyclonal antibody against human LONP2, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Targets peroxisomal Lon protease, useful for studying peroxisomal protein degradation and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	LONP2 Rabbit pAb
Remarks/Alias	LONP, PLON, LONPL, PSLON, LONP2
Species	Human
Gene ID	83752
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human LONP2 (NP_113678.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86WA8
Protein Weight	95kDa
Shipping	Ice bag