

Anti-Loricrin Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P23490

Item Number: CM0023232

Introduction

Loricrin rabbit polyclonal antibody validated for WB, IF/ICC, IF-P, ELISA; reacts with human, mouse, rat. Targets major keratinocyte envelope protein. Ideal for dermatology and keratinization studies.

[Learn More](#)

Parameter	Value
Product Name	LORICRIN Rabbit pAb
Remarks/Alias	LOR; LORICRIN
Species	Human
Gene ID (Human)	4014
Gene ID	4014
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LORICRIN (NP_000418.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23490
Protein Weight	26kDa
Shipping	Ice bag