

Anti-Lox Polyclonal Antibody For Wb - P28300

Item Number: CM0029275

Introduction

Explore this LOX rabbit polyclonal antibody (P28300), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat LOX; recombinant protein immunogen; 47kDa. Ideal for studying collagen/elastin crosslinking and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	LOX Rabbit pAb
Remarks/Alias	AAT10; LOX
Species	Human
GeneID (Human)	4015
GeneID	4015
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human LOX (NP_002308.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28300
Protein Weight	47kDa
Shipping	Ice bag