

Anti-Loxl1 Polyclonal Antibody For Wb, If-P, Elisa - Q08397

Item Number: CM0014866

Introduction

Rabbit polyclonal antibody against human LOXL1, validated for WB, IF-P, and ELISA. Recognizes human and mouse LOXL1, a lysyl oxidase critical for collagen and elastin cross-linking in extracellular matrix.

[Learn More](#)

Parameter	Detail
Product Name	LOXL1 Rabbit Polyclonal Antibody
Remarks/Alias	LOL, LOXL, LOXL1
Species	Human
GeneID (Human)	4016
GeneID	4016
Immunogen	Recombinant fusion protein containing amino acids 95-370 of human LOXL1 (NP_005567.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08397
Protein Weight	63kDa
Shipping	Ice bag