

Anti-Loxl1 Rabbit Monoclonal Antibody For Wb, Elisa - Q08397

Item Number: CM0004547

Introduction

Rabbit monoclonal antibody targeting human LOXL1, suitable for WB and ELISA applications. Validated for cross-reactivity with human samples. Protein weight 63kDa.

[Learn More](#)

Parameter	Value
Product Name	LOXL1 Rabbit mAb
Remarks/Alias	LOL; LOXL; LOXL1
Species	Human
GeneID (Human)	4016
GeneID	4016
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 95-370 of human LOXL1. (NP_005567.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q08397
Protein Weight	63kDa
Shipping	Ice bag